

## The Automated Sample Preparation and Analysis of a Pharmaceutical Cream using the Focus Sample Processing Robot

Diane Nicholas

### Introduction

The extraction of active ingredients from pharmaceutical products, especially creams, can often be a long and tedious process involving many steps and using strongly acidic and basic reagents.

An automated extraction method has been written closely following the customers existing manual method, the main difference being that the volumes have been reduced to allow the use of 10 mL headspace vials. The method uses the Focus Sample Processing Robot for the extraction of an active ingredient (AI) for the quality control of a pharmaceutical cream.

### Instrumentation

- ATASGL Focus Robotic Sample Processor
- Agilent 6890 GC with split/splitless injector and FID
- Agilent ChemStation

The Focus is equipped with the dilutor-dispenser option, which enables the addition of larger volumes of reagents to the sample vial, than the 20  $\mu$ L side-port syringe can hold. The Focus is software controlled and uses magnets to transport the vials between the sample tray and the agitator. Two vials are used per sample, the first contains the sample to which reagents are added to extract the active ingredients. The extracts are then transferred to a second vial and are back-extracted into a solvent that is suitable for injection into the GC. After all extractions are complete the syringe kit is changed and a low volume split injection is made.

### Method

The sample preparation, using the Focus, is a multi-step method. The Focus is loaded with the four reagent vials (alkaline methanol, internal standard (IS), chloroform and acid); two wash solvents (methanol and chloroform); methanol to flush the dilutor; sample vials containing approximately 375 mg of cream; and an equivalent number of second sample vials which are empty.

The Focus adds alkaline methanol solution (2.5 mL) and IS (20  $\mu$ L), transports the vial to the agitator and shakes the sample while it picks-up the chloroform (2.5 mL). After addition of the chloroform the active ingredient is extracted into the methanol and the layer is removed to a second vial. A second extraction of the chloroform is made with alkaline methanol (1.25 mL) and the two extracts are combined. The

Focus transports the second vial, containing the extracts, to the agitator and shakes them before adding firstly acid (500  $\mu$ L) then chloroform (2.5 mL). After all extractions have been performed and the syringe has been changed a 2  $\mu$ L injection is made from the lower chloroform layer. The figure below outlines the process:




Figure: Overview of the pharmaceutical cream extraction method

## *Performance Data*

---

**Table 1:** System suitability calculations

	<b>IS area</b>	<b>AI area</b>	<b>Ratio AI/IS</b>
<b>Cal 1</b>	547.192	1406.862	2.5711
<b>Cal 2</b>	552.111	1413.317	2.5598
<b>Cal 3</b>	558.245	1427.342	2.5568
<b>Cal 4</b>	557.803	1421.862	2.5490
<b>Cal 5</b>	547.438	1410.582	2.5767
<b>Cal 6</b>	554.590	1425.088	2.5696
<b>Average</b>	552.8966	1417.5089	2.5638
<b>Std Dev</b>	4.8684	8.3903	0.0103
<b>%RSD</b>	<b>0.881</b>	<b>0.592</b>	<b>0.403</b>

## *Results*

---

See Table 2 on page 3 for the results of 10 extracted samples and 1 extracted standard.

## *Conclusions*

---

The preliminary results shown here indicate that the Focus system is suitable for the extraction of a pharmaceutical cream, giving %RSDs well far below the limit of 2% and assay results well within the limits of 0.9-1.1 mg/g. Further work will be carried out to determine the robustness of this technique.

## Results

**Table 2:** Results for 10 extracted samples and 1 extracted standard

ID	Sample weight (g)	IS peak area	AI peak area	Std area ratio AI/IS	Smp area ratio AI/IS	Assay result (mg/g)	Mean result (mg/g)	%LS
Cal		558.580	1426.096	2.5531				
Cal		556.927	1427.468	2.5631				
Sample 1	0.39025	580.707	1506.023		2.5934	1.015		101
Sample 1	0.39025	575.757	1494.195		2.5952	1.015	<b>1.02</b>	102
Cal		549.641	1416.939	2.5779				
Cal		552.887	1428.785	2.5842				
Sample 2	0.40164	578.057	1523.065		2.6348	1.002		100
Sample 2	0.40164	573.584	1521.406		2.6525	1.008	<b>1.01</b>	101
Cal		553.106	1431.543	2.5882				
Cal		553.108	1429.087	2.5837				
Sample 3	0.38878	563.515	1448.559		2.5706	1.010		101
Sample 3	0.38878	567.769	1459.442		2.5705	1.010	<b>1.01</b>	101
Cal		556.504	1431.586	2.5725				
Cal		534.145	1383.070	2.5893				
Sample 4	0.37085	573.403	1400.106		2.4417	1.005		101
Sample 4	0.37085	572.828	1395.951		2.4369	1.003	<b>1.00</b>	100
Cal		560.123	1430.782	2.5544				
Cal		564.756	1447.565	2.5632				
Sample 5	0.36315	561.400	1340.135		2.3871	1.004		100
Sample 5	0.36315	576.219	1363.045		2.3655	0.995	<b>1.00</b>	99
Cal		569.292	1459.532	2.5638				
Cal		567.088	1454.506	2.5649				
Sample 6	0.35901	573.068	1348.141		2.3525	1.001		100
Sample 6	0.35901	571.738	1347.061		2.3561	1.002	<b>1.00</b>	100
Cal		560.789	1436.828	2.5622				
Cal		555.699	1429.016	2.5716				
Sample 7	0.38715	591.867	1449.859		2.4496	0.966		97
Sample 7	0.38715	593.659	1453.460		2.4483	0.966	<b>0.97</b>	97
Cal		549.803	1417.701	2.5786				
Cal		550.739	1435.541	2.6066				
Sample 8	0.38784	564.631	1446.760		2.5623	1.009		101
Sample 8	0.38784	572.087	1457.320		2.5474	1.003	<b>1.01</b>	100
Cal		561.169	1444.675	2.5744				
Cal		552.708	1425.642	2.5794				
Sample 9	0.37168	600.103	1458.046		2.4297	0.998		100
Sample 9	0.37168	583.384	1418.301		2.4312	0.999	<b>1.00</b>	100
Cal		561.931	1442.681	2.5674				
Cal		565.071	1452.572	2.5706				
Sample 10	0.39156	584.051	1506.032		2.5786	1.006		101
Sample 10	0.39156	587.847	1515.846		2.5786	1.006	<b>1.01</b>	101
Cal		559.632	1442.107	2.5769				
Cal		602.012	1497.729	2.4879				

**K Factor: 3.495**  
**% RSD (all calibrants): 0.8590**